

Technical Information

MitoFlamma[®] Green

Overview

MitoFlamma[®] Green is a mitochondria-selective green fluorescent dye that allows to detect mitochondrial morphology in living cells. MitoFlamma[®] Green appears to preferentially accumulate in mitochondria regardless of mitochondrial membrane potential in certain cell types, making it a possible tool for determining mitochondrial mass. MitoFlamma[®] Green enables researchers to observe mitochondrial activity, localization and abundance as well as monitoring the effect of drugs or other external stimuli on the mitochondrial function. MitoFlamma[®] Green probes are retained in the mitochondria during the fixation step and after permeabilization with some detergents during subsequent processing steps. After fixation, labeled samples can be applied in a variety of experiments such as immunocytochemistry, *in situ* hybridization, microplate-based analysis, etc.




Figure 1. Excitation/emission spectra of MitoFlamma[®] Green




Figure 2. HeLa cell staining with MitoFlamma[®] Green

Left: DAPI, middle: MitoFlamma Green, right: merged image

Before beginning of mitochondrial labeling Experiment

Materials and equipment required but not provided

- DMSO
- PBS buffer or suitable growth medium for live cell imaging
- Aldehyde based fixatives such as paraformaldehyde for cell fixation
- Aldehyde based detergents such as Triton® X-100
- Micropipette
- Fluorescence microscope
- 37 °C incubator

Experimental protocols

Labeling of mitochondria in fixed or suspended cells with MitoFlamma® Green

1. Prepare the fixed/suspended cells by seeding appropriate cells onto confocal dishes 6 hours before.
2. Adjust the temperature of the media to 37 °C.
3. Prepare a 1 µM labeling MitoFlamma® Green solution by diluting original 1 mM DMSO solution with the media from step 3, *protect from light.
4. Aspirate the culture medium from the cells grown dishes (for suspended cells, centrifuge the dishes in advance), and add 1 µM labeling MitoFlamma® Green solution to the dishes and incubate for 1 h at 37 °C under 5% CO₂ in humid condition.
*the incubation may be varied for optimal labeling condition.
5. Aspirate the labeling media, and rinse the cells twice with PBS at 37 °C.
6. Analyze stained cells with appropriate equipment

Fixation (continuing from step 5 in labeling experiment)

1. Fix the cells for 15 min in 4% paraformaldehyde at 37 °C.
2. Aspirate the fixative, and rinse the cell twice with PBS.
3. View the cells using FITC filter set.

Permeabilization (optional)


- When permeabilization is needed for subsequent steps such as immunocytochemistry, incubate fixed cells in buffer containing detergent such as 0.2% Triton® X-100.
- Following permeabilization, rinse the cells in buffer and proceed with immunocytochemistry procedure.
- Alternatively, the cells may be permeabilized by incubating in ice-cold acetone for 5 minutes, and then washed in PBS.
- Even when cells are not going to be labeled with an antibody, this acetone-permeabilization step may be useful in improving the signal-to background ratio.

Custom Labeling Service


Based on accumulated know-how and technologies, BioActs provide a wide range of custom services such as protein fluorescence labeling, organic synthesis, oligonucleotide synthesis upon customers' request. Our reliable technology has acknowledged by our clients from domestic and overseas universities, institutions, in vitro diagnostic and pharmaceutical companies and has enabled to steadily conduct their requirements. In addition, we can introduce fluorescent materials to many other compounds such as organic and inorganic compounds, drugs, hormones, polymer, peptides, proteins, antibodies, etc. We also can provide chemical and optical analytical data, along with cell and animal experiments.



Nucleic acid



Peptide/Protein



Antibody



Small molecules
/Polymer

Technical Support

ADDRESS

BioActs CO., LTD. DK Tower 10TH F., 595 beon-gil 9, Cheongneung-daero, Namdong-gu, Incheon, 21666, Korea

PHONE & FAX

Tel: +82-32-818-9100

Fax: +82-32-818-8206

WEBSITE

<http://www.bioacts.com/>

MAILS

order@bioacts.com (Order Support)

support@bioacts.com (Customer Support)

ivd@bioacts.com (B2B/Bulk Order Support)

SDS (Safety Data Sheets) You can find SDS at www.bioacts.com, the official website of BioActs.

CoA (Certificate of Analysis) provides detailed quality information of each product. To see CoA, check the lot number written on each product's page at www.bioacts.com, when having trouble with check, contact to our technical support team

Copyright 2009-2017 BioActs All rights reserved. This information is subject to change without notice.

This product can be used only for research purpose, and it can't be used to treat or diagnose human body or animals.

The product cannot be resold to other regions except for South Korea.

We guarantee that this product was manufactured in compliance with spec in this document. This product is for researchers who have received professional education and training, and the guarantee is valid only when the appropriate users used it for correct purpose and method. The guarantee applies to the first purchaser only, and cannot be granted to others. All the models and samples provided to the purchaser are examples of general forms and status, and cannot be consistent. BioActs limits the compensation for inappropriate products to replacement or refund, and this guarantee is not responsible for 1) accidents, disasters or irresistible problems not caused by defect in materials or technology 2) misuse, mistake or carelessness of the user 3) damage due to not observing the instructions regarding product usage 4) damage from mixed use with other defective products 5) products that come from unofficial distribution path. BioActs does not guarantee that the product is free of defect. Other than cases defined by the guarantee, within the maximum range permitted by the proper law, BioActs is not responsible for direct, indirect, special, random or consequential damage that results from violating the guarantee or conditions, or from legal theories. Also, the damages BioActs are not responsible for include but are not limited to damage during use, damage of profit, damages of loss of actual or expected profit, damage regarding money spending, business loss, opportunity loss, credit loss, reputation loss, and indirect or direct damage that occur due to failure of all actions or results used with BioActs products.