

Technical Information

ICG NHS ester labeling kit

Overview

ICG NHS ester Labeling Kit is designed for effective antibody or protein labeling. ICG (Indocyanin green) dye is clinically approved near infrared (NIR) fluorescent dye. NIR fluorescence allows to observe the deep image from the surface of skin and being utilized in a wide range of research fields. The maxima of Ex/Em values are at 785/811nm. ICG might be excited using 750-800 nm laser line or LED and displays excellent optical property. NHS esters readily react with primary amine on antibody or protein by forming stable amide bond and, ipso facto, it is most popular and effective functional group incorporated into fluorescent reagent for protein or antibody labeling. ICG NHS ester labeling kit contains all required contents for antibody and protein labeling. We offer ICG NHS ester labeling kit as a convenient method which enable the labeling process to be easy.

Table 1. Components of ICG NHS ester Labeling Kit

Part No.	Contents	Amount	Storage
ILSM01	Component A, ICG NHS ester	6 μ L	-20 $^{\circ}$ C
ILSM02	Component B, Reaction buffer	500 μ L	2-8 $^{\circ}$ C
ILSM03	Component C, Wash buffer	2.5 mL	2-30 $^{\circ}$ C
ILSM04	Component D, Filtration tube	3 tubes	




Figure 1. NHS ester reaction scheme for chemical conjugation to a primary amine.

Stable amide bond is formed and N- hydroxysuccinimide is released after the reaction.

Application

ICG NHS ester labeling kit enable detection of biomarkers. ICG-conjugated antibody can specifically visualize the target cells expressing the receptor. For example, peptides, proteins, and antibodies are cancer-specific ligands that form biological conjugation with ICG through ICG NHS ester labeling kit protocol. And ICG labeling kit can be used for labeling of FDA-approved monoclonal antibodies such as Panitubumab, Daclizumab, or Trastuzumab. The ICG-conjugated monoclonal antibodies are switch optical imaging probes

that activated only in target cells.

Before beginning of labeling Experiment

Materials and equipment required but not provided

- Micropipette
- 37 °C incubator
- Microcentrifuge
- Microtubes
- PBS or appropriate buffer

Experimental protocols

Labeling of protein or antibody with ICG NHS ester labeling kit

1. Add 100 µL of wash buffer and sample solution containing 50-200 µg protein.
2. Centrifuge at 8,000 rcf for 15 min.
3. Add 100 µL of reaction solution and 1.74 µL of ICG NHS ester.
4. Mix by pipetting and incubate at 37°C for 10 min
5. Add 100 µL of wash buffer
6. Centrifuge at 8,000 rcf for 15 min.
7. Add 200 µL of wash buffer.
8. Centrifuge at 8,000 rcf for 15 min.
9. Repeat step 7 and 8. A_{800}
10. Resuspend with 200 µL of PBS or appropriate buffer.
11. Transfer the resultant to microtube and store at 4°C.

Filtered solution in the collection tube after centrifugation can be discarded.

Calculation of dye / protein ratio

- Ratio (dye / protein) =
$$\frac{A_{800}/147,000}{(A_{280}-A_{800}\times 0.075)/(\epsilon \text{ of protein})}$$

A800: absorbance at 800 nm

A280: absorbance at 280 nm

ϵ : molar absorption coefficient of protein at 280 nm

- ϵ of IgG is 216,000
- ϵ of ICG in PBS is 147,000

Custom Labeling Service


Based on accumulated know-how and technologies, BioActs provide a wide range of custom services such as protein fluorescence labeling, organic synthesis, oligonucleotide synthesis upon customers' request. Our reliable technology has acknowledged by our clients from domestic and overseas universities, institutions, in vitro diagnostic and pharmaceutical companies and has enabled to steadily conduct their requirements. In addition, we can introduce fluorescent materials to many other compounds such as organic and inorganic compounds, drugs, hormones, polymer, peptides, proteins, antibodies, etc. We also can provide chemical and optical analytical data, along with cell and animal experiments.



Nucleic acid



Peptide/Protein



Antibody



Small molecules
/Polymer

Technical Support

ADDRESS

BioActs CO., LTD. DK Tower 10TH F., 595 beon-gil 9, Cheongneung-daero, Namdong-gu, Incheon, 21666, Korea

PHONE & FAX

Tel: +82-32-818-9100

Fax: +82-32-818-8206

WEBSITE

<http://www.bioacts.com/>

MAILS

order@bioacts.com (Order Support)

support@bioacts.com (Customer Support)

ivd@bioacts.com (B2B/Bulk Order Support)

SDS (Safety Data Sheets) You can find SDS at www.bioacts.com, the official website of BioActs.

CoA (Certificate of Analysis) provides detailed quality information of each product. To see CoA, check the lot number written on each product's page at www.bioacts.com, when having trouble with check, contact to our technical support team

Copyright 2009-2017 BioActs All rights reserved. This information is subject to change without notice.

This product can be used only for research purpose, and it can't be used to treat or diagnose human body or animals.

The product cannot be resold to other regions except for South Korea.

We guarantee that this product was manufactured in compliance with spec in this document. This product is for researchers who have received professional education and training, and the guarantee is valid only when the appropriate users used it for correct purpose and method. The guarantee applies to the first purchaser only, and cannot be granted to others. All the models and samples provided to the purchaser are examples of general forms and status, and cannot be consistent. BioActs limits the compensation for inappropriate products to replacement or refund, and this guarantee is not responsible for 1) accidents, disasters or irresistible problems not caused by defect in materials or technology 2) misuse, mistake or carelessness of the user 3) damage due to not observing the instructions regarding product usage 4) damage from mixed use with other defective products 5) products that come from unofficial distribution path. BioActs does not guarantee that the product is free of defect. Other than cases defined by the guarantee, within the maximum range permitted by the proper law, BioActs is not responsible for direct, indirect, special, random or consequential damage that results from violating the guarantee or conditions, or from legal theories. Also, the damages BioActs are not responsible for include but are not limited to damage during use, damage of profit, damages of loss of actual or expected profit, damage regarding money spending, business loss, opportunity loss, credit loss, reputation loss, and indirect or direct damage that occur due to failure of all actions or results used with BioActs products.